

Celda de Oxígeno marca Nuova modelo E-16/0, compatible con:



- ✓ IT Gambert: M-16
- ✓ Maxtec: MAX-16, MAX-41
- ✓ MSA (Catalyst Research): 710644
- ✓ Puritan Bennett 740, 760, 840, PB2000 Test: 4-072214-0
- ✓ Viasys (Carefusion): AVEA/Infant flow CPAP
- ✓ GE: Versamed



Contacto:


E-mail: ventas@inavtek.com

Celular: 55 6686 9108

Celda de Oxígeno marca Nuova modelo E-17/A2, compatible con:



- ✓ Airshields
C2000/C2000e/C200 (Hill-ROM): 8362030, MU13223
- ✓ Atom: 2000, 2100 Incubator, 2100G
- ✓ City Technologies: MOX-01
- ✓ Criticare: 644 oxygen cell, 8100 Poet Plus, 8500 Poet, 8500 Poet IQ, IQ2
- ✓ Datascope: 0600-00-0002
- ✓ Envitec: OOM102, OOM102-C, OOM103
- ✓ Marquette (GE): R-22
- ✓ Hamilton: Oxygen Cell HM-12, Rapahel, Galileo, Arabella, Multniex
- ✓ Heinen und Loewenstein: OOM102
- ✓ Henleys: 5556
- ✓ Hudson: 5566, V-12
- ✓ IT-Gambert: M-12
- ✓ Maxtec: Max-12A
- ✓ Ohmeda (GE) Centiva S/5: 6051-0000-216
- ✓ Sensor TEC (BPR Medical): ST-12
- ✓ Spacelabs: 015-0132-03, 90514
- ✓ Taema: Cesar
- ✓ Teledyne: R22MED, R36,



Contacto:

E-mail: ventas@inavtek.com

Celular: 55 6686 9108

Celda de Oxígeno marca Nuova modelo E-17/J, compatible con:



- ✓ Airshields (Dräger): 6735142
- ✓ Analytical Industries: PSR-11-917-j
- ✓ BCI: R-17
- ✓ Bird: 6400, Vent/8400, Vent/VIP, Vent/Sentry, 66030
- ✓ Breas (GE): Vivo 50
- ✓ Dameca: Alphamed
- ✓ Datex-Engstrom: Erica Elvira
- ✓ Envitec: OOM202-1, OOM103-1
- ✓ F. Stephan: Christina 1107600009 Vent.
- ✓ Hamilton: HM-11
- ✓ Imagen y Medicina: Matisse
- ✓ Invivo: 9445 sensor, Precess 3160,
- ✓ Hudson: 5567, V-11
- ✓ Imagen y Medicina: Matisse
- ✓ Invivo: 9445 sensor, Precess 3160, Millenia 3500, 3155A
- ✓ IT-Gambert: M03
- ✓ Kontron: ABT5100
- ✓ Lifecare: PLV-100/102/37-079
- ✓ Maxtec: Max-13 (CAG13)
- ✓ Medec Belgium: M-03
- ✓ Ohmeda (GE): 6051-0000-217
- ✓ Oxitron: 3310
- ✓ Purittan Bennett 7820, 7200, 850: 4-020933-00
- ✓ Respironics: Esprit / lifecare 1001454
- ✓ Sensidyne: SV11A
- ✓ Sensor TEC (BPR Medical): ST-11
- ✓ SIARE: Morpheus, Perseo
- ✓ Taema: Eole, Cesar
- ✓ Vascular Technology: 100580
- ✓ Ventrex: V-11/A
- ✓ Viasys (Carefusion): 6400 / 8400 / Sensor 60030

Contacto:

E-mail: ventas@inavtek.com

Celular: 55 6686 9108

Celda de Oxígeno marca Nuova modelo E-15/2, compatible con:



- ✓ Analytical Instruments: PSR-11-15/2
- ✓ Ceramatec: CAG-11
- ✓ DP Medical: MVA210
- ✓ Dräger:
- ✓ 6850645/6803290 Cicero, Cato, Babylog, Sulla, Oxycom, Oxydig, Savina, Evita 2 & 4, Evita Dura, Evita XL, PM8, Fabius, Fabius GS, Anemone, Narcomed 1A, 2, 2B, 3, 4, 6850645 oxygen cell
- ✓ Envitec: OOM201
- ✓ Hamilton: HM-13, HM-14
- ✓ Hewlett Packard 15201A
- ✓ Hudson: 5557, V-14
- ✓ IT Gambert: M-02/L2
- ✓ Maxtec: Max-11
- ✓ MSA (Catalyst Research): 655264/655263
- ✓ Ohmeda (GE): 6051-0000-224
- ✓ Oxitron: 3710
- ✓ Paragon Services: PS-2000
- ✓ PPG Biomedical, Irisa: 505055-001
- ✓ Sensor Tec (BPR Medical): ST-14
- ✓ Teledyne: C-61339 (R23MED)
- ✓ Vascular Technology: 101372
- ✓ Ventrex: V-14A

Contacto:

E-mail: ventas@inavtek.com

Celular: 55 6686 9108

Celda de Oxígeno marca Nuova modelo E-75/4, compatible con:



- ✓ Maxtec: Max-250K
- ✓ Viasys (Carefusion): Vela

Contacto:

E-mail: ventas@inavtek.com

Celular: 55 6686 9108

Celda de Oxígeno marca Nuova modelo E-15/4, compatible con:



- ✓ Analytical Industries: PSR-11-95/4
- ✓ Ceramatec: CAG-10
- ✓ IT-Gambert: M-10
- ✓ Maxtec: MAX-10
- ✓ Ohmeda (GE):
- ✓ Smart Vent7900, 7100 Ventilator, Modulus SE, S/5 Aespire, Aestiva, 3000/5, 6050-004-110
- ✓ Sensidyne SV-20A
- ✓ Teledyne R30MED
- ✓ Ventrex V-20A




Contacto:

E-mail: ventas@inavtek.com

Celular: 55 6686 9108

Celda de Oxígeno marca Nuova modelo E-43/0, compatible con:



- ✓ Envitec: OOM112
- ✓ IT-Gambert: M43
- ✓ Ohmeda (GE): Giraffe
6600-1278-600, 6600-0120-850

Contacto:

E-mail: ventas@inavtek.com

Celular: 55 6686 9108